

COLOR CODES OF SELECTED GLOVE SIZES

Gloves

Gloves must be marked with the size indication based on the size of the hand they are to cover, or with the coded color band appropriate to the glove size indicated in prospectus A.1. The prospectus A.2 provides the traditional glove sizes that correspond to the color of gloves encoded in prospectus A.1. Gloves that do not match the dimensions of the prospectus A.1 must not be equipped with bands with these colors. Appendix B contains information on hand measurements and glove sizes that fit different hand sizes.

Prospectus A.1

Size of color-coded gloves (Size in mm)

Band Colour	Size (\pm tolerance in mm)												
	l_0 ± 5	l_1 ± 3	l_2 ± 3	l_3 ± 3	l_4 ± 3	l_5 ± 3	l_6 ± 3	l_7 ± 3	l_8 ± 3	l_9 ± 5	l_{10} ± 5	l_{11} ± 5	l_{12} ± 3
Brown	188	110	46	95	69	79	71	55	8	35	95	40	10
Green	200	114	48	100	72	82	74	58	8	35	95	40	10
White	212	118	50	105	75	85	77	61	9	40	100	45	10
Red	225	122	52	110	78	88	80	64	10	45	105	50	10
Blue	237	127	55	118	81	91	83	67	11	50	110	55	10
Orange	250	132	58	125	84	94	86	70	12	55	115	60	10

Prospectus A.2

Color codes matching traditional glove sizes

<u>Band Colour</u>	<u>Glove Size EU</u>
Brown	5-5.5
Green	6-6.5
White	7-7.5
Red	8-8.5
Blue	9-9.5
Orange	10

Flat Gloves

Six sizes of flat type gloves have been designated by contrast-band bands distinguished by different colors. The A.1 perspective shows the value and tolerance for each dimension of each glove size. The size of the gloves is measured with the help of finger gauges. The glove finger gauges have parallel sides and a rounded end as shown in Figure A.2 so that they can be inserted into the glove finger. They consist of sheets of 1 mm thickness with widths corresponding to the values from b, to above as indicated in prospectus A.3. The gauges are used by introducing them gently into the corresponding finger and measuring the length of the gauge inserted in relation to the corresponding indications given in Figure A.1.

Prospectus A.3

Dimensions of gauges used to measure the size of color-coded gloves

Dimensions in mm

Band Colour	Calipers for Sizes				
	b_1	b_2	b_3	b_4	b_5
Brown	28	27	27	25	22
Green	30	29	29	27	24
White	32	31	31	29	26
Red	35	33	33	31	28
Blue	38	35	36	33	30
Orange	41	37	38	35	32

Figure A.1

Flat glove dimensions from l_0 to l_{12} indicated in prospectus A.1, measured by means of a scale using the calipers from b_1 to b_5 as indicated in prospectus A.3 in relation to gloves marked with color-coded wrist bands indicating the size.

